

DIGITAL BLOOD PRESSURE

FEATURE CHART III



Digital Blood Pressure Summary Chart



ADC Item #	6012	6015	6016	6021	6022	6023
Inflation	Semi-Auto	Automatic	Advanced Auto	Automatic	Advanced Auto	Advanced Auto
Deflation	Auto Valve	Elec. Auto	Elec. Auto	Elec. Auto	Elec. Auto	Elec. Auto
Type	Upper Arm	Wrist	Wrist	Upper Arm	Upper Arm	Upper Arm
Smart Logic	No	Yes	Yes	Yes	Yes	Yes
Pulse	Yes	Yes	Yes	Yes	Yes	Yes
Date/Time	No	Yes	Yes	Yes	Yes	Yes
Memory	Last Reading	30	2 X 30	99	2 X 99	2 X 99
Memory Zone	1	1	2	1	2	2
Average Mode	No	No	Yes (MAM)	No	Yes (MAM)	Yes (MAM)
Sizes Available	1	1	1	3	1	2
Limb Range	9"-13" Adult	5"-7 5/8" Wrist	5"-7 5/8" Wrist	7.5"-9" Sm. Adult 9"-13" Adult 13"-17" Lg. Adult	9"-13" Adult	9"-13" Adult 13"-17" Lg. Adult
Quick Start Button	Yes	Yes	Yes	Yes	Yes	Yes
Cuff Storage Compartment	No	N/A	N/A	Yes	Yes	No
Carry Case	No	Yes	Yes	Yes	Yes	Yes
Battery Type	4 "AA"	2 "AAA"	2 "AAA"	4 "AA"	4 "AA"	4 "AA"
A/C Adapter Available	No	No	No	Yes	Yes	Yes
Latex Content	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free	Latex-free
Warranty	5 Year	5 Year	5 Year	5 Year	5 Year	5 Year
Case Pack	6	24	12	6	6	6
Carton Pack	3	12	3	3	3	3
Catalog Page	36	36	36	38	38	37
BHS Validation	Yes	No	Yes	Yes	Yes	Yes
Replaces Old Unit	6006N	6011	6011	6013	6014	6017
PC Connectivity	N/A	N/A	Yes	N/A	Yes	Yes

Digital Blood Pressure Cross Reference Chart



ADC	6012	6015	6016 ^A	6021	6022 ^{A/3}	6023 ^{A/B}
A & D	UA704	UB328	UB512	N/A	UA774	UA787EJ
A & D	N/A	N/A	UB521	UA853AC	UA767PC	UA851
A & D	N/A	N/A	N/A	UA767	N/A	UA789AC
A & D	N/A	N/A	N/A	UA767 Plus	N/A	N/A
A & D	N/A	N/A	N/A	N/A	N/A	N/A
Lumiscope	1103	1143	1145	1133	1137	N/A
Lumiscope	N/A	N/A	N/A	N/A	N/A	N/A
Mabis	N/A	04-235-001	04-248-001	N/A	04-342-000	04-346-000 ^B
Mabis	04-263-001	N/A	04-255-001	N/A	04-320-001	N/A
Mabis	04-340-000	N/A	04-790-001	04-330-001 ²	04-310-001	N/A
Mabis	N/A	N/A	04-237-001	N/A	N/A	N/A
Mabis	N/A	N/A	N/A	N/A	N/A	N/A
Omron	N/A	N/A	N/A	N/A	N/A	HEM-790IT
Omron	N/A	HEM-629	HEM-650	HEM-711AC ¹	HEM-780	N/A
Omron	N/A	N/A	HEM-670IT	HEM-712C	N/A	N/A
Omron	N/A	N/A	N/A	HEM-711DLX	N/A	N/A
Prestige	HM-10	HM-40	N/A	HM-20	N/A	N/A
Samsung	N/A	N/A	BWD-6007	N/A	N/A	N/A
Samsung	BVM-1007	HD-5007	N/A	N/A	BVA-2007	N/A
Samsung	N/A	N/A	N/A	BSD-3007	N/A	BSP-4007

- * **Note 1** Includes AC Adapter, ADC adapter sold separately
- * **Note 2** High speed model. ADC's all work at "regular" speed
- * **Note 3** Includes ADC PC Software sold separately

Note A Only ADC offers Average mode (averages 3 readings)
Note B ADC unit comes with PC Software

Digital Blood Pressure Model #6012 ADvantage™ Semi-Automatic Blood Pressure Monitor

Features:	Benefits:
Automatic Off	Conserves battery life
Extra Large Display	For easy reading
Single Large Button	Simplifies use
Latex-Free	Reduces risk of allergic reactions
Functional Check	Checks all display functions & calibration self-test
Short Tone Indicates Starting	Patient is aware of measurement
Long Tone Indicates Finished	Patient is aware of completion
Auto-Valve	Simplifies measurement
Digital Inflation Instructions	Guides patient with cuff pressure
Digital Error Messages	Guides patient when there's a problem
Weighs Less Than 2lbs	Lightweight for travel
D-Ring Contoured Cuff	Easy, self-application
Last Reading Memory	Simplifies record keeping
8 Bit Microprocessor	Improved logic and accuracy
5 Year Warranty	Longer consumer protection
Includes Record Book	For convenient record keeping
BHS Validated	*Customer's Assurance of Accuracy*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6012
SEMI-AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6021 ADvantage™ Automatic Blood Pressure Monitor

Features:	Benefits:
Automatic Off	Conserves battery life
Smart Logic	No need to set systolic inflation level
Extra Large Display	For easy reading
Memory Stores Last 99 Readings	Simplifies record keeping
Date/Time Stamp	Keeps more relevant records
Latex-Free	Reduces risk of allergic reactions
Irregular Heartbeat Detection	Signals appearance of irregular heartbeat
Functional Check	Displays battery warning
5 Beeps Indicate Starting	Patient is aware of measurement
Long Tone Indicates Finished	Patient is aware of completion
Digital Error Messages	Guides patient when there's a problem
Available in Small and Large Adult	For small & large adult patients
Cuff Storage Compartment	Keeps unit compact for neatness
A/C Adapter Option	Conserves battery life
D-Ring Contoured Cuff	Easy, self-application
8 Bit Microprocessor	Improved logic and accuracy
5 Year Warranty	Longer consumer protection
Includes Record Book	For convenient record keeping
Deluxe Carry Case	For convenient travel and storage
BHS Validated	*Customer's Assurance of Accuracy*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6021
AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6022
ADvantage™ PLUS Advanced Automatic Blood Pressure Monitor

Features:	Benefits:
Automatic Off	Conserves battery life
Smart Logic	No need to set systolic inflation level
Two Zone Memory	Ideal for use with two patients
Date/Time Stamp	Keeps more relevant records
Latex-Free	Reduces risk of allergic reactions
Irregular Heartbeat Detection	Signals appearance of irregular heartbeat
Functional Check	Displays battery warning
5 Beeps Indicate Starting	Patient is aware of measurement
Long Tone Indicates Finished	Patient is aware of completion
Digital Error Messages	Guides patient when there's a problem
Cuff Storage Compartment	Keeps unit compact for neatness
D-Ring Contoured Cuff	Easy, self-application
A/C Adapter Option	Conserves battery life
Extra Large Display	For easy reading
2 Zone 99 Reading Memory	Simplifies record keeping
PC Software available	For easy data analysis
Average Mode	Improved accuracy
Wide Range Soft Cuff	For medium to large adult patients
8 Bit Microprocessor	Improved logic and accuracy
5 Year Warranty	Longer consumer protection
Includes Record Book	For convenient record keeping
Deluxe Carry Case	For convenient travel and storage
BHS Validated	*Customer's Assurance of Accuracy*

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.



ADVANTAGE™ 6022 PLUS
AUTOMATIC BP MONITOR

Digital Blood Pressure Model #6023 ADvantage™ **ULTRA** Advanced Blood Pressure Monitor

Features:

Benefits:

Automatic Off	Conserves battery life
Smart Logic	Automatically determines ideal inflation levels
Traffic Light Indicator	Compares reading with National Standards
Date/Time Stamp	Keeps more relevant records
Latex-Free	Reduces risk of allergic reactions
Irregular Heartbeat Detection	Signals appearance of irregular heartbeat
Digital Error Messages	Guides patient when there's a problem
Wide Range Rigid Cuff	Easy, self-application for medium and large adults
Advanced 2 Zone Memory	Recalls last 99 readings for two people
MAM Advantage™ Average Mode	Improved accuracy
PC Link Software	Connects monitor to your PC to download and print data, graphs, charts, and tables from readings
5 Year Warranty	Longer consumer protection
AC Adapter	Conserves battery life
USB Cable	For easy connection to your PC
Deluxe Carry Case	Convenient for travel and storage
*BHS Validated	Customer's Assurance of Accuracy

***BHS (British Hypertension Society)** validation was issued to only a handful of manufacturers including ours.

It is equivalent to an endorsement by the AMA.

ADVANTAGE™ 6023 **ULTRA**
AUTOMATIC BP MONITOR



Digital Blood Pressure Model #6015 ADvantage™ Basic Wrist Blood Pressure Monitor



Features:

Automatic Off
Smart Logic
Memory Stores Last 30 Readings with Date/Time Stamp
Latex-Free
Functional Check
Digital Error Messages
Weighs Less Than 1lb
Includes Carry Case
5 Year Warranty
Includes Record Book
Single Button Operation
Measures at Wrist
Contoured Cuff

Benefits:

Conserves battery life
 No need to set systolic inflation level
 Simplifies record keeping
 Reduces risk of allergic reactions
 Displays battery warning
 Guides patient when there's a problem
 Compact & lightweight for travel
 Convenient for travel
 Longer consumer protection
 For convenient record keeping
 Simplifies use
 Improved portability & convenience
 Easy, self-application

Digital Blood Pressure Model #6016 ADvantage™ Advanced Wrist Blood Pressure Monitor



Features:

Automatic Off
Smart Logic
Two Zone Memory
Date/Time Stamp
Latex-Free
Long Tone Indicates Finished
Digital Error Messages
Contoured Cuff
Memory Stores Last 30 Readings Each
PC Software Available
Average Mode
8 Bit Microprocessor
5 Year Warranty
Includes Record Book
Weighs Less Than 1lb
Includes Carry Case
Measures at Wrist
BHS Validated

Benefits:

Conserves battery life
 No need to set systolic inflation level
 Ideal for use with two patients
 Keeps more relevant records
 Reduces risk of allergic reactions
 Patient is aware of completion
 Guides patient when there's a problem
 Easy, self-application
 Simplifies record keeping
 For easy data analysis
 Improved accuracy
 Improved logic and accuracy
 Longer consumer protection
 For convenient record keeping
 Compact & lightweight for travel
 Convenient for travel
 Improved portability & convenience
 Customer's Assurance of Accuracy

BHS (British Hypertension Society) validation was issued to only a handful of manufacturers including ours. Equivalent endorsement to the AMA.



ADC
 55 Commerce Drive
 Hauppauge, NY 11788

tel: 631-273-9600
 toll free: 1-800-232-2670
 fax: 631-273-9659

www.adctoday.com
info@adctoday.com



ADC (UK) Ltd.
 Unit 6, PO14 1TH
 United Kingdom



6 3 4 7 8 2 1 0 5 4 4 7 2
 digbpc2004 - rev 3