

Our latest must-have for busy professionals:

# ADTEMP™ MINI 432

## Non-Contact Thermometer


ADC's non-contact Adtemp Mini 432 is compact, lightweight, and fast. Featuring an ingenious design and incredible portability, the Adtemp 432 is sure to become an essential instrument for busy professionals.

### NO SKIN CONTACT REQUIRED!

Minimize the risk of cross contamination



Retractable, **breakaway lanyard** included – the perfect instrument for patient rounds!



- One-second readings with 9-reading memory
- Compact, lightweight size: measures 3.5" and weighs 1.6 oz
- Color-coded LCD display with fever warning
- °F/°C switchable display
- Up to 1,000 measurements on a single charge
- Good for 300,000 measurements, lifetime
- Micro USB cable included
- Clips to wearable, retractable, breakaway lanyard, included
- Audible beep when measurement is complete
- 2-year warranty

### COLOR-CODED FEVER WARNING



Normal



Elevated



Fever

55 Commerce Drive  
Hauppauge, New York 11788  
Tel: 1-800-ADC-2670  
Fax: 631-273-9659  
[www.adctoday.com](http://www.adctoday.com)

**NOW AVAILABLE!**

[www.adctoday.com/products/432](http://www.adctoday.com/products/432)

